\\USER\FMRIF\[XT-ID:93-M-0170]Renzo\20250219_SRA_templ\scout_sag

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further	On
preparation	
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.25 ms
TE	1.53 ms
Averages	1
Concatenations	1
Filter	B1 filter
Coil elements	AC

Contrast - Common

TR	3.25 ms
TE	1.53 ms
Flip angle	16.0 deg

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution - Common

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

Resolution - iPAT

Reference scan mode	Integrated
·	

Resolution - Filter Image

Image Filter	Off	
Distortion Corr.	Off	
Prescan Normalize	Off	
Normalize	Off	
B1 filter	On	
Unfiltered images	Off	

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.25 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
Р	0.0 mm
Н	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	Н
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off

System - Miscellaneous

Coil Select Mode Default	
--------------------------	--

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Flip angle	16.0 deg
Measurements	1
Time to center	6.3 s

Inline - Inline

Subtract	Off	
Measurements	1	
StdDev	Off	
Save original images	On	

Inline - MIP

MIP-Sag	Off	
MIP-Cor	Off	
MIP-Tra	Off	
MIP-Time	Off	
Save original images	On	

Inline - Composing

Inline Composing	Off
Distortion Corr	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	16.0 deg
Measurements	1
Contrasts	1

Inline - MapIt

TR	3.25 ms
TE	1.53 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

Mode Off

\\USER\FMRIF\[XT-ID:93-M-0170]Renzo\20250219_SRA_templ\scout_trans

TA: 0:20 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further	On
preparation	
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	22.2 %
Slices per slab	144
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.25 ms
TE	1.53 ms
Averages	1
Concatenations	1
Filter	B1 filter
Coil elements	AC

Contrast - Common

TR	3.25 ms
TE	1.53 ms
Flip angle	16.0 deg

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution - Common

FoV read	260 mm	
FoV phase	100.0 %	
Slice thickness	1.6 mm	
Base resolution	160	
Phase resolution	100 %	
Slice resolution	69 %	
Phase partial Fourier	6/8	
Slice partial Fourier	6/8	
Trajectory	Cartesian	

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

Resolution - iPAT

Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	On
Unfiltered images	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice oversampling	22.2 %
Slices per slab	144
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.25 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 A30.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Initial Position	L0.0 A30.0 H0.0
L	0.0 mm
Α	30.0 mm
Н	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	Н
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off

System - Miscellaneous

Coil Select Mode Default	
--------------------------	--

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm	
Excitation	Non-sel.	

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Flip angle	16.0 deg
Measurements	1
Time to center	8.3 s

Inline - Inline

Subtract	Off	
Measurements	1	
StdDev	Off	
Save original images	On	

Inline - MIP

MIP-Sag	Off	
MIP-Cor	Off	
MIP-Tra	Off	
MIP-Time	Off	
Save original images	On	

Inline - Composing

Inline Composing	Off	
Distortion Corr	Off	

Inline - MapIt

Save original images	On
Maplt	None
Flip angle	16.0 deg
Measurements	1
Contrasts	1

Inline - MapIt

TR	3.25 ms
TE	1.53 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

Mode	Off	

$\verb|\USER\FMRIF|[XT-ID:93-M-0170]| Renzo | 20250219_SRA_templ | VASO_fMRI_AP_test|| | 20250219_SRA_test|| | 202$

TA: 0:47 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further	Off
preparation	
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P40.6 H21.7 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1
•	

Contrast - Dynamic

,	
Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s

Resolution - Common

FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
Base resolution	226
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

	_
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P40.6 H21.7 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group 1 Position R7.2 P40.6 H21.7 mm Orientation T > C-44.4 Phase enc. dir. A >> P AutoAlign Initial Position R7.2 P40.6 H21.7 R 7.2 mm P 40.6 mm H 21.7 mm Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4 > S 0.0		
Orientation T > C-44.4 Phase enc. dir. A >> P AutoAlign Initial Position R7.2 P40.6 H21.7 R 7.2 mm P 40.6 mm H 21.7 mm Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4	Slab group	1
Phase enc. dir. A >> P AutoAlign Initial Position R7.2 P40.6 H21.7 R 7.2 mm P 40.6 mm H 21.7 mm Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4	Position	R7.2 P40.6 H21.7 mm
AutoAlign Initial Position R7.2 P40.6 H21.7 R 7.2 mm P 40.6 mm H 21.7 mm Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4	Orientation	T > C-44.4
Initial Position	Phase enc. dir.	A >> P
R 7.2 mm P 40.6 mm H 21.7 mm Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4	AutoAlign	
P 40.6 mm H 21.7 mm Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4	Initial Position	R7.2 P40.6 H21.7
H 21.7 mm Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4	R	7.2 mm
Initial Rotation 0.00 deg Initial Orientation T > C T > C -44.4	P	40.6 mm
Initial Orientation $T > C$ T > C -44.4	Н	21.7 mm
T > C -44.4	Initial Rotation	0.00 deg
	Initial Orientation	T > C
> S 0.0	T > C	-44.4
	> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	Н
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

PATRef FA RF duration RF BWT product RF BWT product Ernst T1 1200 ms PATRef prep. shots 10 Volume dummy shots 0 Dummy Measurements 0 ETL per RTEB 1 Invert PE Off Min. TE if PF Con Echo Time Shift On NORDIC SVDPC On Sym VASO Off Invert 3D Off Disable PF reco Disable PF reco Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs FIDN		
RF BWT product Ernst T1 1200 ms PATRef prep. shots 10 Volume dummy shots 0 Dummy Measurements 0 ETL per RTEB 1 Invert PE Off Min. TE if PF On Echo Time Shift On Ramp Sampling On NORDIC On SVDPC On Sym VASO Off Invert RO Off Invert RO Off Diasble PF reco Off Save sampling Off PE VComp Off Water Excnone- External PC per Series Saturation RF per Shot FIDNavs -none- EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On GRAPPA Regularization 10 10^-6 HSN RF power scale Inversion 10 0 ms Relaxation Delay 500 ms	PATRef FA	3 deg
Ernst T1 1200 ms PATRef prep. shots 10 Volume dummy shots 0 Dummy Measurements 0 ETL per RTEB 1 Invert PE Off Min. TE if PF On Echo Time Shift On Ramp Sampling On NORDIC On SVDPC On Sym VASO Off Invert RO Off Invert RO Off Invert RO Off Invert RO Off Save sampling Off Save sampling Off Vater Exc. External PC per Series Saturation RF per Shot FIDNavs -none-EPI rise time factor Modify Ice Config GRAPPA Regularization Hold Policy Inversion Delay Relaxation Delay Relaxation Delay Some Automatical Policy Inversion Delay Relaxation Delay Relaxation Delay Some Interest Into Off One Disable PF reco Off One Disable PF reco Off Off Off Off Off Off Off Off Off Of	RF duration	2500 us
PATRef prep. shots Volume dummy shots Dummy Measurements ETL per RTEB Invert PE Min. TE if PF Echo Time Shift Ramp Sampling NORDIC SVDPC Sym VASO Dual-pol. EPI Invert RO Invert 3D Disable PF reco Disable PF reco Save sampling Off Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Modify Ice Config GRAPPA Regularization HSN RF power scale Inversion Inversion Delay Relaxation Delay Relaxation Delay Relaxation Delay Off Double Description On On On On On On On On On	RF BWT product	8
Volume dummy shots Dummy Measurements Dummy Measure	Ernst T1	1200 ms
Dummy Measurements ETL per RTEB Invert PE Min. TE if PF Cho Time Shift Ramp Sampling NORDIC SVDPC On Sym VASO Off Dual-pol. EPI Invert RO Invert 3D Off Disable PF reco Off Save sampling Off Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay Relaxation Delay Off Off Off Off Off Off Off O	PATRef prep. shots	10
ETL per RTEB Invert PE Off Min. TE if PF Cho Time Shift On Ramp Sampling On NORDIC SVDPC On Sym VASO Off Dual-pol. EPI Invert RO Invert 3D Off Disable PF reco Disable PF reco Off Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Mosaic DICOMs Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay Relaxation Delay On	Volume dummy shots	0
Invert PE Min. TE if PF Echo Time Shift On Ramp Sampling On NORDIC SVDPC Sym VASO Off Dual-pol. EPI Invert RO Invert 3D Off Disable PF reco Off Save sampling Off Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Mosaic DICOMs Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay On On On On On On On On On O	Dummy Measurements	0
Min. TE if PF Echo Time Shift Ramp Sampling On NORDIC SVDPC On Sym VASO Off Dual-pol. EPI Invert RO Invert 3D Off Disable PF reco Disable PF reco Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Mosaic DICOMs Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay Ron On	ETL per RTEB	1
Echo Time Shift Ramp Sampling On NORDIC SVDPC SVDPC On Sym VASO Off Dual-pol. EPI Invert RO Invert 3D Off Disable PF reco Disable PF reco Off Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Mosaic DICOMs Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay On SVDPC On	Invert PE	Off
Ramp Sampling NORDIC SVDPC Sym VASO Off Dual-pol. EPI Invert RO Invert 3D Disable PF reco Disable PF reco Off Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Mosaic DICOMs Modify Ice Config GRAPPA Regularization HSN RF power scale Inversion Delay Relaxation Delay On On On On On On On On On O	Min. TE if PF	On
NORDIC SVDPC SVM VASO Off Dual-pol. EPI Invert RO Invert 3D Disable PF reco Disable PF reco Off Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs EPI rise time factor Mosaic DICOMs Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay Off Off On Off On	Echo Time Shift	On
SVDPC Sym VASO Off Dual-pol. EPI Invert RO Invert 3D Off Disable PF reco Off Disable PF reco Off Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs FIDNavs EPI rise time factor Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay Off Off Off Off Off Off Off Off Off Of	Ramp Sampling	On
Sym VASO Off Dual-pol. EPI Invert RO Off Invert 3D Off Disable PF reco Off Disable PF reco Off Save sampling Off PE VComp Off Water Excnone- External PC per Series Saturation RF per Shot FIDNavs -none- EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	NORDIC	On
Dual-pol. EPI Invert RO Invert 3D Off Invert 3D Off Disable PF reco Off Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs FIDNavs FIDNavs EPI rise time factor Modify Ice Config G-factor map GRAPPA Regularization Inversion Delay Relaxation Delay Off Off Off Off Off Off Off Off Off Of	SVDPC	On
Invert RO Invert 3D Off Disable PF reco Off Disable PF reco Off Save sampling Off PE VComp Off Water Exc. External PC Saturation RF FIDNavs FIDNavs FIDNavs FIDNavs FIDNavs FIDNavs On Modify Ice Config GRAPPA Regularization Inversion Delay Relaxation Delay Off Off Off Off Off Off Off Off Off Of	Sym VASO	Off
Invert 3D Off Disable PF reco Off Disable PF reco Off Save sampling Off PE VComp Off Water Excnone- External PC per Series Saturation RF per Shot FIDNavs -none- EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	Dual-pol. EPI	Off
Disable PF reco Disable PF reco Off Save sampling Off PE VComp Off Water Exc. External PC Saturation RF FIDNavs FIDNavs FIDNavs FIDNavs FIDNavs On Modify Ice Config G-factor map GRAPPA Regularization Inversion Delay Relaxation Delay Off Off Off Off Off Off Off Off Off Of	Invert RO	Off
Disable PF reco Save sampling Off PE VComp Water Exc. External PC Saturation RF FIDNavs FIDNavs FIDNavs FIDNavs On Modify Ice Config GRAPPA Regularization HSN RF power scale Inversion Delay Relaxation Delay Off Off Off Off Off Off Off Off Off Of	Invert 3D	Off
Save sampling PE VComp Off Water Exc. External PC Saturation RF FIDNavs FIDNavs FIDNavs FIDNavs On Modify Ice Config GRAPPA Regularization HSN RF power scale Inversion Delay Relaxation Delay Off Off Off Off Off Off Off Off Off Of	Disable PF reco	Off
PE VComp Water Excnone- External PC per Series Saturation RF per Shot FIDNavs -none- EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^-6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	Disable PF reco	Off
Water Excnone- External PC per Series Saturation RF per Shot FIDNavs -none- EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^-6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	Save sampling	Off
External PC Saturation RF FIDNavs FIDNavs -none- EPI rise time factor Mosaic DICOMs Modify Ice Config G-factor map GRAPPA Regularization HSN RF power scale Inversion Delay Relaxation Delay Residue Shot Per Series per Series none- Per Series Per Series None- Per	PE VComp	Off
Saturation RF per Shot FIDNavs -none- EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^-6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	Water Exc.	-none-
FIDNavs -none- EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^-6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	External PC	per Series
EPI rise time factor 1.10 Mosaic DICOMs On Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^-6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	Saturation RF	per Shot
Mosaic DICOMs Modify Ice Config G-factor map GRAPPA Regularization HSN RF power scale Inversion Delay Relaxation Delay On On On On On On On On On O	FIDNavs	-none-
Modify Ice Config On G-factor map Off GRAPPA Regularization 10 10^-6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	EPI rise time factor	1.10
G-factor map GRAPPA Regularization HSN RF power scale Inversion Delay Relaxation Delay Off 10 10^-6 3.00 500 ms Relaxation Delay 0 ms	Mosaic DICOMs	On
GRAPPA Regularization 10 10^-6 HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	Modify Ice Config	On
HSN RF power scale 3.00 Inversion Delay 500 ms Relaxation Delay 0 ms	G-factor map	Off
Inversion Delay 500 ms Relaxation Delay 0 ms	GRAPPA Regularization	10 10^-6
Relaxation Delay 0 ms	HSN RF power scale	3.00
,	Inversion Delay	500 ms
Var. FA /MAGEC 4	Relaxation Delay	0 ms
	Var. FA /MAGEC	4

Mode	Off

TA: 12:35 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further	Off
preparation	
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1
·	·

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	202
Pause after meas.	0.0 s

Resolution - Common

FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
Base resolution	226
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off

Resolution - Common

Interpolation Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Initial Position	R7.2 P41.2 H17.8
R	7.2 mm
P	41.2 mm
Н	17.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-44.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Set-n-Go Protocol	Off	
Table position	Н	

Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10^-6
HSN RF power scale	3.00
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Mode	Off

$\verb|\USER\FMRIF|[XT-ID:93-M-0170]Renzo|| 20250219_SRA_templ|| VASO_fMRI_AP_run2|| 20250219_SRA_templ|| 20250219_SR$

TA: 12:35 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	202
Pause after meas.	0.0 s

Resolution - Common

FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
Base resolution	226
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off

Resolution - Common

Interpolation Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Initial Position	R7.2 P41.2 H17.8
R	7.2 mm
P	41.2 mm
Н	17.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-44.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Set-n-Go Protocol	Off
Table position	Н

Table position	0 mm	
Inline Composing	Off	

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10^-6
HSN RF power scale	3.00
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Mode	Off

$\verb|\USER\FMRIF|[XT-ID:93-M-0170]Renzo|| 20250219_SRA_templ|| VASO_fMRI_AP_run3|| 20250219_SRA_templ|| VASO_fMRI_AP_run3|| 20250219_SRA_templ|| 20250219_SRA$

TA: 12:35 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	202
Pause after meas.	0.0 s

Resolution - Common

FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
Base resolution	226
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off

Resolution - Common

Interpolation Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Initial Position	R7.2 P41.2 H17.8
R	7.2 mm
P	41.2 mm
Н	17.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-44.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Set-n-Go Protocol	Off
Table position	Н

Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

ocquerioe - opeciai	
PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
SVDPC	On
Sym VASO	Off
p	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
. =	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
	1.10
	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10^-6
HSN RF power scale	3.00
	500 ms
	0 ms
Var. FA /MAGEC	4

Mode	Off

$\verb|\USER\FMRIF\[XT-ID:93-M-0170]| Renzo | 20250219_SRA_templ | VASO_fMRI_PA| | A superiordization | A superiordiz$

TA: 0:47 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	P >> A
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

Contrast - Dynamic

,	
Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s

Resolution - Common

FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
Base resolution	226
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	P >> A
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	P >> A
AutoAlign	
Initial Position	R7.2 P41.2 H17.8
R	7.2 mm
P	41.2 mm
Н	17.8 mm
Initial Rotation	180.00 deg
Initial Orientation	T > C
T > C	-44.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off	
Table position	Н	
Table position	0 mm	
Inline Composing	Off	

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10^-6
HSN RF power scale	3.00
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Mode	Off	
IVIOUC	Oli	

$\verb|\USER\FMRIF|[XT-ID:93-M-0170]Renzo|| 20250219_SRA_templ\\| T1234_AP|| T123$

TA: 3:38 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5l

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further	Off
preparation	
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L0.0 P6.5 F3.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	
Slices per slab	186
FoV read	190 mm
FoV phase	108.6 %
Slice thickness	0.82 mm
TR 1	28.5 ms
TR 2	52560 ms
TE 1	8.27 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	28.5 ms
TR 2	52560 ms
TE 1	8.27 ms
Multi-echo spacing	19 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1093.5 ms
TI 2	2860.5 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	14

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	4
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s

Resolution - Common

FoV read	190 mm	
FoV phase	108.6 %	
Slice thickness	0.82 mm	
Base resolution	232	
Phase resolution	100 %	
Slice resolution	100 %	
Phase partial Fourier	6/8	

Resolution - Common

Slice partial Fourier	Off	
Interpolation	Off	

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	48
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L0.0 P6.5 F3.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slices per slab	186
FoV read	190 mm
FoV phase	108.6 %
Slice thickness	0.82 mm
TR 1	28.5 ms
TR 2	52560 ms

Geometry - AutoAlign

Slab group	1
Position	L0.0 P6.5 F3.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	
Initial Position	L0.0 P6.5 F3.3
L	0.0 mm
P	6.5 mm
F	3.3 mm
Initial Rotation	44.40 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Set-n-Go Protocol	Off
Table position	Н
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.29 ms
Bandwidth	938 Hz/Px

Sequence - Part 2

EPI factor	14
Segmentation	14
RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Turbo factor	62

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
Ernst T1	1200 ms

Sequence - Special

PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	On
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Mode	Off
------	-----

$\verb|\USER\FMRIF|[XT-ID:93-M-0170]Renzo|| 20250219_SRA_templ|| VASO_fMRI_AP_run4|| AP_run4|| AP_r$

TA: 12:35 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further	Off
preparation	
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	202
Pause after meas.	0.0 s

Resolution - Common

FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
Base resolution	226
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off

Resolution - Common

Interpolation Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Initial Position	R7.2 P41.2 H17.8
R	7.2 mm
P	41.2 mm
Н	17.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-44.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Set-n-Go Protocol	Off
Table position	Н

Table position	0 mm	
Inline Composing	Off	

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10^-6
HSN RF power scale	3.00
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Mode	Off

$\verb|\USER\FMRIF|[XT-ID:93-M-0170]Renzo|| 20250219_SRA_templ\\| T1234_PA|| T123$

TA: 3:38 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5l

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L0.0 P6.5 F3.3 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	
Slices per slab	186
FoV read	190 mm
FoV phase	108.6 %
Slice thickness	0.82 mm
TR 1	28.5 ms
TR 2	52560 ms
TE 1	8.27 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	28.5 ms
TR 2	52560 ms
TE 1	8.27 ms
Multi-echo spacing	19 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1093.5 ms
TI 2	2860.5 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	14

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	4
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s

Resolution - Common

FoV read	190 mm
FoV phase	108.6 %
Slice thickness	0.82 mm
Base resolution	232
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

Resolution - Common

Slice partial Fourier	Off	
Interpolation	Off	

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	48
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	L0.0 P6.5 F3.3 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
Slices per slab	186
FoV read	190 mm
FoV phase	108.6 %
Slice thickness	0.82 mm
TR 1	28.5 ms
TR 2	52560 ms

Geometry - AutoAlign

Slab group	1
Position	L0.0 P6.5 F3.3 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	
Initial Position	L0.0 P6.5 F3.3
L	0.0 mm
P	6.5 mm
F	3.3 mm
Initial Rotation	-135.60 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Set-n-Go Protocol	Off
Table position	Н
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	Н
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.29 ms
Bandwidth	938 Hz/Px

Sequence - Part 2

EPI factor	14
Segmentation	14
RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Turbo factor	62

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
Ernst T1	1200 ms

Sequence - Special

PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	On
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Mode Off	
----------	--

TA: 12:35 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1
·	·

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	202
Pause after meas.	0.0 s

Resolution - Common

_		
Ī	FoV read	180 mm
	FoV phase	99.1 %
	Slice thickness	0.80 mm
	Base resolution	226
	Phase resolution	100 %
	Slice resolution	100 %
	Phase partial Fourier	6/8
l	Slice partial Fourier	Off

Resolution - Common

Interpolation Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Initial Position	R7.2 P41.2 H17.8
R	7.2 mm
P	41.2 mm
Н	17.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-44.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Set-n-Go Protocol	Off
Table position	Н

Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10^-6
HSN RF power scale	3.00
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

		_
Mode	Off	

\\USER\FMRIF\[XT-ID:93-M-0170]Renzo\20250219_SRA_tempI\VASO_fMRI_AP_run6

TA: 12:35 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	56.9 ms
TR 2	3688 ms
TE 1	20.00 ms
Multi-echo spacing	48 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1192.8 ms
TI 2	2558.4 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	202
Pause after meas.	0.0 s

Resolution - Common

FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
Base resolution	226
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off

Resolution - Common

Interpolation Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	80
Acc. factor 3D	2
Ref. lines 3D	24
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	4

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off	
Elliptical filter	Off	

Geometry - Common

Slab group	1
Slabs	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	180 mm
FoV phase	99.1 %
Slice thickness	0.80 mm
TR 1	56.9 ms
TR 2	3688 ms

Geometry - AutoAlign

Slab group	1
Position	R7.2 P41.2 H17.8 mm
Orientation	T > C-44.4
Phase enc. dir.	A >> P
AutoAlign	
Initial Position	R7.2 P41.2 H17.8
R	7.2 mm
P	41.2 mm
Н	17.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-44.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	Fat sat.

Set-n-Go Protocol	Off
Table position	Н

Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	Н
Table position	0 mm
MSMA	S-C-T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R9.8 P37.5 H22.1 mm
! Orientation	T > C-41.9
! Rotation	90.00 deg
! R >> L	195 mm
! A >> P	175 mm
! F >> H	38 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.160541 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.12 ms
Bandwidth	1006 Hz/Px

Sequence - Part 2

EPI factor	42
Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

Sequence - Special

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10^-6
HSN RF power scale	3.00
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Mode	Off